



HOSPITECH™
HEALTHCARE FOR HUMANITY



RENAL CARE

HOSPITECH MANUFACTURING SERVICES SDN BHD
199601034155 (406507-H)

AUREUMAEX MEDICARE MARKETING (M) SDN BHD
198501011707 (144159-P)

COMPANY HISTORY

Hospitech is an integrated medical devices provider specialising in the design, development, value-added manufacturing and distribution of medical devices. In addition, Hospitech is also an independent service provider of Ethylene Oxide (EtO) sterilization, primarily for medical devices.

Hospitech Marketing Sdn Bhd (HMKT) was established in 1985 as a distribution company of generic medical devices and critical care devices. During the early years, HMKT started only with a small range of medical products, focused primarily on generic medical devices. In anticipation of an increase in local demands for generic medical disposables device, Hospitech Manufacturing Services Sdn Bhd (HMFG) was formed and became its manufacturing division in 1998. Subsequently, in 2000, HMFG commenced manufacturing of critical care medical devices.

In 2009 change of ownership by Private Equity which composed of Proventeus which hold 65% of share. Currently, 100% of share the ownership structure have been hold by Proventeus and there are 348 employees employed under medical division companies. On 2014, Hospitech Marketing Sdn Bhd (HMKT) have change its name to Aureumaex Medicare Marketing (M) Sdn Bhd (AMM).

Presently, Hospitech manufactures a large family of disposable medical devices which consist of generic and critical care medical devices under the HOSPITECH™ brand. Hospitech has an EtO sterilisation facility which is ISO 11135:2014 certified. This facility is supported by its fully automated in-house developed EtO Control System (ECS).

Hospitech is currently operating from Sub Lot 1, Lot 6536, Batu 6, Jalan Bukit Kemuning, 40460 Shah Alam, Selangor Darul Ehsan, Malaysia. This 2.2 acres of premise houses Hospitech's local and international marketing team, EtO plant, warehouse as well as sales and administrative office.

Vision

To be the No 1 and most reputable healthcare manufacturer in Malaysia with a highly effective organization that offer great place to work and provide self-esteem, commitment and fulfillment.

Mission

Commitment to excel quality, speed and cost through innovation and continues improvement.

Quality Policy

Consistently provide high quality healthcare products and services.

Meet or exceed customer's expectation in term of safety, performance and cost.

Implementation of continuous improvement program and maintain effectiveness of Quality Management System.

QUALITY CERTIFICATES

Hospitech has obtained its certifications for ISO 13485: 2016 from SGS United Kingdom Ltd Systems Services Certification (SGS UK) as well as CE mark for some of its medical devices.



RENAL CARE

content

- Hemodialysis Blood Tubing Set (DEHP Free) Pg 4
- Dialyzer (DEHP Free) Pg 6
- Double/Triple Lumen Pg 8
- AV Fistula Needle Pg 9
- Hemodialysis Bicarbonate (Catridge) Pg 10
- Citric Acid Pg 10
- Priming Bag Pg 11
- Surgical Spirit Pg 11
- Sodium Hypochlorite Pg 11






HEMODIALYSIS BLOOD TUBING SET (DEHP FREE)

Description:

- Model for all hemodialysis machine and for all configurations.
- Single peel with ETO Sterilisation.
- Ready for use. For single use only.
- Each chamber have a monitor line which is color coded.
- The arterial and venous bloodline have a blood sampling port which is color coded.
- The sampling port is made of latex free material.
- The tubing material shall be made from medical grade PVC (DEHP free), supporting document provided.
- The configuration may subjected to customer requirement.

Features:

1. Pump segment for conventional and high flow rates
 - Flanges keep the pump segment in the correct position.
2. One-stop blood clamp
 - Avoids the hazards of multistop clamp.
3. Priming / Rinseback connector
 - Priming set can be transferred directly to the patient connector simplifying rinseback.
4. Unique luer lock patient connector
 - Twist-off features removes even the tightest fitting fistula needle set. Affixing ISO standard.
5. Easy to handle tubing
 - Soft, remains kink-free. Frosted surface prevents sticking.
6. Compact one piece filter
 - Primes quickly. Low profile allows the chamber to be placed high in the air detector, hereby reducing false alarm.
7. Finger grips
 - Secure connectors are your patients lifeline. Only Hospitech's connectors have finger grips for easier connects and disconnects - even with wet gloves. Affixing ISO standard.
8. Priming recirculation connector
 - Vented for priming and with luer lock recirculation.
9. Self sealing injection sites
 - Extremely soft injection site material assures leak free injection / aspiration even when using very large or poorly siliconized needles.
10. Models for all hemodialysis machines and for all configurations
 - Drip chambers are available to meet all requirement
 - Innovative design reduces foaming and clotting at either high or low flow rates.



DIALYZER (DEHP FREE)

PPS

Description:

- Unique process technology influences separation profile of membrane formation.
- Thicker separation layer and sharpened sieving profile.
- Optimized dialysate side distribution.
- Higher mechanical strength.
- Resulting low albumin loss and high molecule removal.
- Maximum clearance of low molecular.
- Very high CIS retention.
- Single pack, ready to use.



SHM-DLPES-P1.8MF-PPS

| PRODUCT CODE | DESCRIPTION | PACKAGING |
|----------------------|---------------------------------------|-------------------|
| PPS MID FLUX | | ALL CTN OF 24 PCS |
| SHM-DLPES-P1.4MF-PPS | DIALYZER PES PLUS 1.4 MIDDLE FLUX-PPS | |
| SHM-DLPES-P1.6MF-PPS | DIALYZER PES PLUS 1.6 MIDDLE FLUX-PPS | |
| SHM-DLPES-P1.8MF-PPS | DIALYZER PES PLUS 1.8 MIDDLE FLUX-PPS | |
| SHM-DLPES-P2.0MF-PPS | DIALYZER PES PLUS 2.0 MIDDLE FLUX-PPS | |
| PPS HIGH FLUX | | |
| SHM-DLPES-P1.4HF-PPS | DIALYZER PES PLUS 1.4 HIGH FLUX-PPS | |
| SHM-DLPES-P1.6HF-PPS | DIALYZER PES PLUS 1.6 HIGH FLUX-PPS | |
| SHM-DLPES-P1.8HF-PPS | DIALYZER PES PLUS 1.8 HIGH FLUX-PPS | |
| SHM-DLPES-P2.0HF-PPS | DIALYZER PES PLUS 2.0 HIGH FLUX-PPS | |
| SHM-DLPES-P2.2HF-PPS | DIALYZER PES PLUS 2.2 HIGH FLUX-PPS | |



SHM-DLPES-P1.8HF-PPS



SHM-DLPES-P1.6MF-PPS



PH

Description:

- Competitive to other synthetic low flux-membrance in terms of hemocompatibility and performance.
- Superior to conventional low-flux membranes in terms of middle molecular removal.
- Higher UFR compared to other semisynthetic and synthetic low flux membranes.
- Single pack, ready to use.

| PRODUCT CODE | DESCRIPTION | PACKAGING |
|---------------------|--------------------------------------|-------------------|
| PH LOW FLUX | | ALL CTN OF 24 PCS |
| SHM-DLPES-P1.8LF-PH | DIALYZER PES PLUS 1.8 LOW FLUX - PH | |
| PH HIGH FLUX | | |
| SHM-DLPES-P1.6HF-PH | DIALYZER PES PLUS 1.6 HIGH FLUX - PH | |
| SHM-DLPES-P1.8HF-PH | DIALYZER PES PLUS 1.8 HIGH FLUX - PH | |
| SHM-DLPES-P1.8HF-PH | DIALYZER PES PLUS 2.0 HIGH FLUX - PH | |

DOUBLE/TRIPLE LUMEN

Description:

- A double lumen catheter has 2 lumens and triple lumen catheter has 3 lumens.
- Allow us to run several different infusion with only one access site.
- Often inserted through a short central line with a wider inner diameter.
- Single peel pack, sterile, ready for use.
- For single use only.




| PRODUCT CODE | DESCRIPTION | | PACKAGING |
|---------------------|---|-----------|---------------|
| DOUBLE LUMEN | | | |
| SL-DSCVC-12X24 | DOUBLE LUMEN STRAIGHT 12FR X 24CM | STRAIGHT | CTN OF 25 PCS |
| SL-DSCVC-12X20 | DOUBLE LUMEN STRAIGHT 12FR X 20CM | STRAIGHT | CTN OF 25 PCS |
| SL-DSCVC-12X15 | DOUBLE LUMEN STRAIGHT 12FR X 15CM | STRAIGHT | CTN OF 25 PCS |
| SL-DCCVC-12X20 | DOUBLE LUMEN CURVE 12FR X 20CM | CURVE | CTN OF 25 PCS |
| SL-DCCVC-12.5X15 | DOUBLE LUMEN CURVE SZ: 12.5FR X 15CM | CURVE | CTN OF 50 PCS |
| SL-DPCCVC-11.5X12 | DOUBLE LUMEN PRE CURVE 11.5FR X 12CM | PRE CURVE | CTN OF 40 PCS |
| TRIPLE LUMEN | | | |
| SL-TLHDC-1215 | TRIPLE LUMEN CURVE CATH 12FRX15CM | CURVE | CTN OF 25 PCS |
| SL-TLHDS-1215 | TRIPLE LUMEN HD STRAIGHT CATH 12FRX15CM | STRAIGHT | CTN OF 25 PCS |
| SL-TLHDS-1220 | TRIPLE LUMEN HD STRAIGHT CATH 12FRX20CM | STRAIGHT | CTN OF 25 PCS |

AV FISTULA NEEDLE

Description:

- Clamp-snap-type allowing simple termination of blood flow.
- Tube-selected elastic material offer the most suitable flexibility to fit its function.
- Wing - Flexible wing to ensure easy turning for the appropriate direction after puncture procedure and for maximum blood flow.
- Color coded wings distinguish gauge size easily and clearly.
- Needle - Ultra thin-wall needle ensures sufficient blood flow. Most advanced siliconized coating to prevent blood coagulation and hemolysis.
- Back Eye - The introduction of an oval vent offer minimum tissue damage and maximum blood flow.
- Tube Length : 300mm. Promotes healthcare worker safety through the prevention of accidental needle stick injury.
- Single peel pack, Sterile, Ready for use.



A.V. Fistula Safety Needle



SHM-AVF-16G



SHM-AVF-15G

| PRODUCT CODE | DESCRIPTION | PACKAGING |
|-------------------|-------------------------------------|----------------------|
| SHM-AVF-15G | A.V. FISTULA NEEDLE SZ:15G (HMS) | BOX OF 50 PCS |
| SHM-AVF-16G | A.V. FISTULA NEEDLE SZ:16G (HMS) | BOX OF 50 PCS |
| SHM-AVF-17G | A.V. FISTULA NEEDLE SZ: 17G (HMS) | BOX OF 50 PCS |
| SHM-AVFS-15GX25MM | A.V. FISTULA SAFETY NEEDLE 15GX25MM | BOX OF 50 PCS |
| SHM-AVFS-16GX25MM | A.V. FISTULA SAFETY NEEDLE 16GX25MM | BOX OF 50 PCS |
| SHM-AVFS-17GX25MM | A.V. FISTULA SAFETY NEEDLE 17GX25MM | BOX OF 50 PCS |
| SHM-AVF-15G SM | ARTERIAL VENOUS FISTULA | 50'S/Box - 500'S/Ctn |
| SHM-AVF-16G SM | ARTERIAL VENOUS FISTULA | 50'S/Box - 500'S/Ctn |
| SHM-AVFR-15G | AVF NEEDLE ROTATING 15G | 50'S/Box - 500'S/Ctn |
| SHM-AVFR-16G | AVF NEEDLE ROTATING 16G | 50'S/Box - 500'S/Ctn |
| SHM-AVFR-17G | AVF NEEDLE ROTATING 17G | 50'S/Box - 500'S/Ctn |
| SHM-AVFR15-15G | AVF NEEDLE ROTATING 15G (15CM) | 50'S/Box - 500'S/Ctn |
| SHM-AVFR15-16G | AVF NEEDLE ROTATING 16G (15CM) | 50'S/Box - 500'S/Ctn |

HEMODIALYSIS BICARBONATE CATRIDGE

Description:

- Helps the hemodialysis equipment clean the blood and balance acidity.
- Used as a dialysate buffer in hemodialysis treatment.
- Comes with weight of 650g.
- Cartridge packaging.
- Single use.



| PRODUCT CODE | DESCRIPTION | PACKAGING |
|--------------|----------------------|-----------|
| HB001 | BICARB CATRIDGE 650G | 10's/CTN |

CITRIC ACID

Description:

- Molecular formula: $C_6H_8O_7$
- Appearance: White, Solid, Crystalline
- Solubility: Soluble in alcohol
- Oder: Odorless
- Specific Gravity: 1.090 - 1.120
- Boiling Point: Decomposes before boiling
- pH(value) : 1.7(100g/1)

Uses:

- Citric acid is used for the cleaning of dialysis equipment.



| PRODUCT CODE | DESCRIPTION | PACKAGING |
|--------------|-------------------------|-----------------------|
| SR-012 | CITRIC ACID (ANHYDROUS) | 1KG/BTL , CTN OF 9 KG |

PRIMING BAG

Description:

- Priming Bag is used in the Haemodialysis bloodline tubing set.
- For single use only.



| PRODUCT CODE | DESCRIPTION | PACKAGING |
|----------------|------------------------------|------------------------|
| SHM-PB-2000ML | PRIMING BAG 2L - NON STERILE | 100'S/PKT, 5'PKTS/CTN. |
| SHM-PBS-2000ML | PRIMING BAG 2L - STERILE | 100'S/CTN |

SURGICAL SPIRIT

Features:

- Used for cleansing the skin before injections or before surgery.
- To use, apply with a clean cotton ball or gauze to clean intended surface.
- Safe and effective for intended use.
- For External Use Only.



| PRODUCT CODE | DESCRIPTION | PACKAGING |
|--------------|--------------------------------|----------------|
| SR-003 | SURGICAL SPIRIT (70%) 4 LITERS | 4 LITRE/BOTTLE |

SODIUM HYPOCHLORITE

Description:

- Liquid Bleach
- Appearance: Clear, light yellow liquid
- Chemical Formula: NACIO
- Solubility: Completely soluble in water
- Physical state at room temperature: Liquid
- Oder: Chlorine-like
- Specific Gravity: 1.090 - 1.120
- Boiling Point: >212F
- pH(value) : >11.80

Packing:

- The material should be packed in air tight plastic containers or as loose tanker.



| PRODUCT CODE | DESCRIPTION | PACKAGING |
|--------------|--------------------------------|----------------|
| SM-011 | SODIUM HYPOCHLORITE 5.25% (5L) | 5 LITRE/BOTTLE |



HOSPITECH™
HEALTHCARE FOR HUMANITY

**Hospitech Manufacturing Services Sdn Bhd,
Aureumaex Medicare Marketing (M) Sdn Bhd,**
Sublot 1, Lot 6536, Batu 6, Jalan Bukit Kemuning,
40460 Shah Alam, Selangor Darul Ehsan,
Malaysia.

Email: enquiry@hospitech.net

Website: www.hospitech.my

Tel: 03-5122 5199

Fax: 03-5122 5204

Kuching Branch

Ground Floor,

SL 6 Of Parent Lot 2855,

Block 7, Muara Tebas Land District,

Tanjong Sebarang Pending Point,

Sejingkat, 93050 Kuching,

Sarawak.

Tel: 082-552 221

More Branches:

Butterworth, Penang | Kota Bharu, Kelantan | Ayer Keroh, Melaka

Johor Bahru, Johor | Batu Pahat, Johor

